



Ref. TS53000S



CE
0459



28/11/2025 19:55:43 – For information purposes only, can be subject to change

Ref. TS53000S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215486	HS Code	90189010
Country	France		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	25
Packed unit / Length (cm)	16	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES