

SARGASSE LONG-SLEEVED WHITE COAT, SIZE 0

Ref. WSSAR07_94



SARGASSE LONG-SLEEVED WHITE COAT, SIZE 0

Ref. WSSAR07_94



Description

- Mixed blouse
-
- Short sleeves
-
- Tailored collar
-
- Twill
-
- 2 low pockets, 1 chest pocket and 1 inner pocket
-
- Closure on the front by snap buttoning
-
- Height clothing +/- 112 cm
-
- 65% polyester, 35% cotton
-
- 200 g / m²
- Available in size 1 (WSSAR07_91), size 2 (WSSAR07_92), size 4 (WSSAR07_94)

SARGASSE LONG-SLEEVED WHITE COAT, SIZE 0

Ref. WSSAR07_94



Features

CE Marking	NONE	Medical device	No
Bar code	3661040203867	HS Code	62114310030
Country	Tunisia		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Packed unit / Height (cm)	1
Packed unit / Length (cm)	37	Packed unit / Weight	0.45
Bare Unit / Weight (kg)	0.45		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Bleaching	Bleaching prohibited	Washing	Wash at 80°C
Model	Mixed model	Cleaning	Dry cleaning
Packed unit / Quantity	1	Size	Size 4