

BLACK NYLON SOFT BAG SPENGLER, MEDIUM

Ref. 535020



Description

Nylon pouch 18x13x6 cm with éclair® zip closure

01/12/2025 21:41:19 – For information purposes only, can be subject to change

BLACK NYLON SOFT BAG SPENGLER, MEDIUM

Ref. 535020



Features

CE Marking	NONE	Medical device	No
Bar code	3700446032395	HS Code	90330000
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---