

GREY SPENGLER (ML) CUFF (22-42CM) FOR TENSORIC & AUTOTENSIO

Ref. 527464



Description

Adult grey cuff (M/L) (22-42cm) for SPENGLER TENSORIC and AUTOTENSIO electronic blood pressure monitor

02/12/2025 22:53:55 – For information purposes only, can be subject to change

GREY SPENGLER (ML) CUFF (22-42CM) FOR TENSORIC & AUTOTENSIO



Ref. 527464

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041908	HS Code	90262040
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.102

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---