

DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSONIC

Ref. 527500P



CE
0459



02/12/2025 21:31:03 – For information purposes only, can be subject to change

DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSONIC



Ref. 527500P

SPENGLER precision in electronic version

Description

Display of 6 Tensonic arm Adult (M/L) Blueberry

-

-VERY HIGH ACCURACY OF
MEASUREMENTS•

SPENGLER measurement• Arrhythmia detector•

Tricolor WHO risk indicator•

3 average calculations available to monitor the average of morning values

EASY TO USE•

Fully automatic

- Large backlit screen with tricolor WHO risk indicator.

- New armband 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL

INFORMATION• Comes with 4 LR6 AA batteries and a carrying case•

Optional AC adapter

- 2 cuff sizes (S, M/L) available as an option

DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC



Ref. 527500P

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446023041	HS Code	90189010
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.59
Bare Unit / Weight (kg)	2.04		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---