

ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC

Ref. 527500



CE
0459

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
3 ANS
3
3 YEARS
GUARANTEE

ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC



Ref. 527500

SPENGLER precision in electronic version

Description

The TENSORIC arm blood pressure monitor is designed for regular therapeutic monitoring. -

-VERY HIGH ACCURACY OF
MEASUREMENTS•

SPENGLER measurement• irregular pulse rate detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values

EASY TO

USE• Fully automatic

• Large backlit screen with tricolor WHO risk indicator.

• New armband 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL

INFORMATION• Comes with 4 LR6 AA batteries and a carrying case•

Optional AC adapter

• 2 cuff sizes (S, M/L) available as an option

ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC



Ref. 527500

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446022624	HS Code	90189010
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.765
Bare Unit / Weight (kg)	0.34		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---