

DISPLAY OF 6 UNITS TENSonic WRIST - WRIST BPM

Ref. 527522P



02/12/2025 21:05:08 – For information purposes only, can be subject to change

DISPLAY OF 6 UNITS TENSORIC WRIST - WRIST BPM



Ref. 527522P

Ideal for regular self-measurement, with a very high level of precision and simple operation

Description

Display of 6 tensiometers on the wrist TENSORIC®.
Availability 2nd quarter for the new color: Blueberry.

VERY HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values

EASY TO

USE• Fully automatic•

Large backlit screen

FAMILY USE•

2 x 60 time-stamped memories

ADDITIONAL INFORMATION

- Wrist circumference 13.5-21 cm
- Comes with 2 LR03 AAA batteries and a storage case

DISPLAY OF 6 UNITS TENSORIC WRIST - WRIST BPM



Ref. 527522P

Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446023058 | HS Code | 90189010 |
| Warranty (months) | 36 | Country | China |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|------|---------------------------|-------|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 1.815 |
| Bare Unit / Weight (kg) | 0.72 | | |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | | | |
|------------------------|---|--|--|
| Packed unit / Quantity | 1 | | |
|------------------------|---|--|--|