


## AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS105

# DDM




CE



02/12/2025 22:30:53 – For information purposes only, can be subject to change

## AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS105



Those lower limb cuffs are autoclavable!

### Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

#### Dimensions of the bladder in cm

Lower limb cuff infant	35.5 x 6 cm
Lower limb cuff child	53 x 6 cm
Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85.5 x 12 cm

**AUTOCLAVABLE POLYESTER LOWER LIMB  
CUFF STRAIGHT 1 BLADDER FOR DDM  
PNEUMATIC TOURNIQUET****Ref. GCS105****Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203674	HS Code	90189084
Warranty (months)	24	Country	France

**DIMENSIONS AND WEIGHT**

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

**SPECIFIC STORAGE CONDITIONS**

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

**OTHER FEATURES**