

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS

Ref. SPG340P-



CE
0459



SPENGLER precision in electronic version

Description

Display of 6 AUTOTENSIO® wrist model. For occasional self-monitoring. Very easy-to-use. - Quick and quiet inflation method - 120 memories with date and time - LCD screen - Fully automatic - Purple cuff - Wrist circumference 13,5 - 21 cm (5,3-8,5 inch). Supplied with 2 LR03 AAA batteries and a storage case.

02/12/2025 22:55:28 – For information purposes only, can be subject to change

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 UNITS



Ref. SPG340P-

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446023027	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.7
Bare Unit / Weight (kg)	1.368		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--