

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL

Ref. SPG345P



CE
0459

DISPOSITIF MÉDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE

AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR DISPLAY OF 6 CORAL



Ref. [SPG345P](#)

SPENGLER precision in electronic version

Description

Display of 6 Autotensio coral wrist.

MEASUREMENT

ACCURACY• SPENGLER MEASUREMENT•

WHO risk indicator

EASY &

CONVENIENT• Fully automatic• 2

buttons• LCD display

• 120 time-stamped memories•

Fast and quiet

ADDITIONAL INFORMATION•

Wrist circumference 13.5-21 cm

• Comes with 2 LR03 AAA batteries and a storage case

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041601	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.62
Bare Unit / Weight (kg)	1.368		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--