

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)

Ref. SPG440P



CE
0459



02/12/2025 22:37:05 – For information purposes only, can be subject to change

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)



Ref. SPG440P

Ease of use in everyday life for self-measurement family use

Description

Display of 6 blood pressure monitors AUTOTENSIO Blueberry.

HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2

buttons•

New cuff 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION

- Comes with 4 LR6 AA batteries and a carrying case• Optional AC adapter
- 2 cuff sizes (S, M/L) available as an option

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM)

Ref. SPG440P



Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	370044602034	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.73
Bare Unit / Weight (kg)	2.148		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--