

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL

Ref. SPG445P



CE
0459



02/12/2025 22:19:39 – For information purposes only, can be subject to change

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL



Ref. SPG445P

Ease of use in everyday life for self-measurement family use

Description

Display of 6 Autotensio arms Coral.

HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2

buttons•

New cuff 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION

- Comes with 4 LR6 AA batteries and a carrying case• Optional AC adapter
- 2 cuff sizes (S, M/L) available as an option

DISPLAY 6 UNITS AUTOTENSIO ARM-ARM BPM ADULT MEDIUM (22X36CM) CORAL

Ref. [SPG445P](#)



Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041588	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.73
Bare Unit / Weight (kg)	2.148		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---