

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL

Ref. SPG445



CE
0459



02/12/2025 22:24:23 – For information purposes only, can be subject to change

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG445

Ease of use in everyday life for self-measurement family use

Description

HIGH MEASUREMENT

ACCURACY•

SPENGLER MEASUREMENT• Arrhythmia detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values EASY TO USE•

Fully automatic• 2

buttons•

New cuff 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL INFORMATION

- Comes with 4 LR6 AA batteries and a carrying case• Optional AC adapter
- 2 cuff sizes (S, M/L) available as an option

AUTOTENSIO UPPER ARM (M/L) BLOOD PRESSURE MONITOR SPENGLER CORAL



Ref. SPG445

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446041571	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	10
Packed unit / Length (cm)	17	Packed unit / Weight	0.594
Bare Unit / Weight (kg)	0.358		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	70
Minimum humidity	5	Maximum humidity	93

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--