

## MONOBLOC CONICAL ADULT ARM CUFF WITH ONE BLADDER AND A BUCKLE FOR DDM TOURNIQUETS

Ref. GBCMS105



The monobloc conical arm cuff with one bladder and a buckle !

### Description

This **Adult L** monobloc arm cuff made of **PU** is decontaminable. With a **conical design** and only **one bladder** it is made c  
a monobloc unit with a **buckle** to secure the attachment.

Width of the bladder : 10cm

Length of the bladder : 47cm

Total length : 70cm

The arm cuff can not be sterilized. It can be dry-cleaned or immersed with paying attention to maintain he connector of the tube out of water.

It is also available as lower limb cuffs

# MONOBLOC CONICAL ADULT ARM CUFF WITH ONE BLADDER AND A BUCKLE FOR DDM TOURNIQUETS



Ref. [GBCMS105](#)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468208129	HS Code	90189084
Country	France		

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.3

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES