

CUFF AND LOWER LIMB CUFF WITH 2 BLADDERS FOR DDM TOURNIQUET

Ref. [GBS202](#)

Cuffs and lower limb cuff with 2 bladders for DDM pneumatic tourniquets !

Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **2 removable silicon bladders** for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff child	2 x (35.5 x 6 cm)
Cuff adult	2 x (59 x 6 cm)
Cuff adult	2 x (69 x 9 cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203513	HS Code	90189084
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES

02/12/2025 21:21:54 – For information purposes only, can be subject to change