

## STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507251



04/12/2025 22:07:41 – For information purposes only, can be subject to change

## STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507251



The Magister II is a high-performance professional stethoscope for versatile auscultation.

### Description

#### HAUTE QUALITÉ ACOUSTIQUE

- Membrane flottante multifréquences
- Pression légère = écoute des basses fréquences (vasculaires et cardiaques)
- Pression ferme = écoute des hautes fréquences (pulmonaires)
- Tubulure simple avec un conduit haute performance
- Lyres en aluminium

#### CONFORT

- Bague antifroid
- Embouts auriculaires en silicone
- Nouveau design de lyres plus adapté à l'anatomie
- Indicateur de pavillon

Disponible en version Silver (507251), Black (507261), Myrtille (507262) et Corail (507252)

#### CONDITIONS DE STOCKAGE :

Pression atmosphérique : 700-1060 hPa

Température : -10° à 40° C

Humidité : 95% max

**STETHOSCOPE MAGISTER II DUAL, SILVER****Ref. 507251****Features**

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 3700446029968 | HS Code                   | 90181990 |
| Warranty (months)           | 60            |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | 16            | Packed unit / Height (cm) | 4        |
| Packed unit / Length (cm)   | 32            | Packed unit / Weight      | 0.43     |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | -10           | Max temperature           | 40       |
| Minimum humidity            | NA            | Maximum humidity          | 95       |
| OTHER FEATURES              |               |                           |          |
| Colors                      | Gray          | Packed unit / Quantity    | 1        |