


VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE

Ref. 518113



CE
0459



04/12/2025 21:39:50 – For information purposes only, can be subject to change

**VAQUEZ-LAUBRY CLASSIC/VELCRO
CUFF/LARGE ADULT (L) / NAVY BLUE****Ref. 518113**

=====

Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. [518113](#)

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008871	HS Code	90189010
Warranty (months)	60	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.802
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		