

ANTIBACTERIAL AND HYDROPHOBIC FILTER FOR ASPIRET ASKIR 20 / 30 EMIVAC

Ref. BCASPI96

CA-MI



Description

- Anti-bacterial and hydrophobic filter
- Diam 65mm / Connector 8mm

04/12/2025 21:38:55 – For information purposes only, can be subject to change

ANTIBACTERIAL AND HYDROPHOBIC FILTER FOR ASPIRET ASKIR 20 / 30 EMIVAC

CA-MI

Ref. BCASPI96

Features

CE Marking	NONE	Medical device	No
Bar code	8054619101209	HS Code	90192000
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.05
Bare Unit / Weight (kg)	0.05		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		