



Ref. MONI011E



CE



04/12/2025 22:12:33 – For information purposes only, can be subject to change

Ref. MONI011E

## Description

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 6948309002582 | HS Code        | 9018193090 |
| Warranty (months) | 6             | Country        | China      |

### DIMENSIONS AND WEIGHT

|                           |     |                           |      |
|---------------------------|-----|---------------------------|------|
| Packed unit / Width (cm)  | nc  | Packed unit / Height (cm) | nc   |
| Packed unit / Length (cm) | nc  | Packed unit / Weight      | 0.15 |
| Bare Unit / Weight (kg)   | 0.1 |                           |      |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 55 |
| Minimum humidity    | 10  | Maximum humidity | 93 |

### OTHER FEATURES