


EAR REUSABLE CLIP SPO2 SENSOR FOR VITALIMATE MINI

Ref. MONI010E



CE
0123



EAR REUSABLE CLIP SPO2 SENSOR FOR VITALIMATE MINI

Ref. [MONI010E](#)

The ear clip for Mtalimate MINI of Spengler

Description

Ear clip SpO2 sensor

To be sued with Vitalimate MINI monitor of Spengler.

- length 1m
- 6 pin sensor

Also available for Vitalimate MINI, are the silicon and rigid reusable sensors in every sizes: adult, peditaric and neonat.

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002568	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES