

Ref. 336402



Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

09/12/2025 20:41:05 – For information purposes only, can be subject to change

Ref. 336402



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028527	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.33
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

09/12/2025 20:41:05 – For information purposes only, can be subject to change