

Ref. 336412

 Spengler



CE



09/12/2025 20:47:05 – For information purposes only, can be subject to change

Ref. 336412



Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028510	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.184

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

09/12/2025 20:47:05 – For information purposes only, can be subject to change