


ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201_R



09/12/2025 20:24:10 – For information purposes only, can be subject to change

ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201_R



Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 20 x 10 x 2 cm

Available colors:

- Azure- Blue

-

Golden

- Grey

- Red

- Green

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067162	HS Code	90184990
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10.2
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.2
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20
Packed unit / Weight	0.074	Bare Unit / Weight (kg)	0.073

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Red	Packed unit / Quantity	1
--------	-----	------------------------	---