



## ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201\_V



09/12/2025 21:05:16 – For information purposes only, can be subject to change

## ALUMINIUM TRAY COLOR 20 X 10 CM, AZUR

Ref. AAX201\_V



### Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 20 x 10 x 2 cm

Available colors:

- Azure- Blue

-

Golden

- Grey

- Red

- Green

### Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3661040067179 | HS Code        | 90184990 |
| Warranty (months) | 12            | Country        | France   |

#### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | 10    | Bare Unit / Width (cm)  | 10    |
| Packed unit / Height (cm) | 2     | Bare Unit / Height (cm) | 2.2   |
| Packed unit / Length (cm) | 20    | Bare Unit / Length (cm) | 20    |
| Packed unit / Weight      | 0.074 | Bare Unit / Weight (kg) | 0.073 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|        |       |                        |   |
|--------|-------|------------------------|---|
| Colors | Green | Packed unit / Quantity | 1 |
|--------|-------|------------------------|---|