



## IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081\_R



## IDEAL+ SINGLE HEAD ADULT, BLUE

Ref. DH7081\_R



Single-head adult stethoscope, designed for reliable and accurate auscultation.

### Description

- Stethoscope with single adult metal horn (Ø: 45 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in double pavilion for adults and a pediatric model

Available in Black (DH7081\_N) or Red (DH7081\_R)

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000137	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1