

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M13000

DDM



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
1 AN
1
YEAR
GUARANTEE



09/12/2025 20:54:14 – For information purposes only, can be subject to change

THE SINGLE USE PRESSURE CUFF EASY FUSE DDM

Ref. M13000



The single use pressure cuff made in France !

Description

The pressure infusion cuffs are medical devices designed for daily use in operating theaters, Intensive Care units or emergency transport vehicles. They can be used to manually accelerate infusion of blood and solutions when gravity alone cannot achieve the desired flow rates.

Easyfuse pressure cuffs are single use and compatible with different sizes of blood/solution bags: 500cc, 1000cc and 3000cc.

Optimized monitoring

- **Single use**, no cleaning and no disinfection
- blood/solution bag clearly visible
- pressure indicator with automatic overpressure regulator
- Controlled decompression facilitated by 3 ways valve
- Résistant as made of TPU

All our pressure infusion cuffs are made in France in our Castelsulier (Agen) manufacturing site..

**THE SINGLE USE PRESSURE CUFF EASY
FUSE DDM**

Ref. M13000

**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203872	HS Code	98109084
Warranty (months)	12	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES