

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03

**ELITE BAGS**



📅 09/12/2025 20:47:46 – For information purposes only, can be subject to change

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03



**Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner**

### Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening
- Material: Polyester gray bitone
- Weight: 0.12 kg.
- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03\_N), blue jeans (MEDIAB03\_B) or pink (MEDIAB03\_RC)
- Sold empty

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03



### Features

|                   |               |                |               |
|-------------------|---------------|----------------|---------------|
| CE Marking        | NONE          | Medical device | No            |
| Bar code          | 8436564542199 | HS Code        | 4202.12.99.90 |
| Warranty (months) | 24            | Country        | China         |

#### DIMENSIONS AND WEIGHT

|                           |       |                         |       |
|---------------------------|-------|-------------------------|-------|
| Packed unit / Width (cm)  | 19    | Bare Unit / Width (cm)  | 10    |
| Packed unit / Height (cm) | 6     | Bare Unit / Height (cm) | 6     |
| Packed unit / Length (cm) | 10    | Bare Unit / Length (cm) | 17    |
| Packed unit / Weight      | 0.248 | Bare Unit / Weight (kg) | 0.012 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|        |      |                        |   |
|--------|------|------------------------|---|
| Colors | Gray | Packed unit / Quantity | 1 |
|--------|------|------------------------|---|