

## CRITICAL URGENCY ELITE BAGS, RED

Ref. MECRIT01\_R

**ELITE BAGS**



11/12/2025 21:35:25 – For information purposes only, can be subject to change

## CRITICAL URGENCY ELITE BAGS, RED

Ref. MECRIT01\_R



### Description

Allows all equipment intended for advanced first aid of emergency services to be transported, facilitating access to equipment quickly and efficiently.

- Interior totally adaptable to the needs of the user thanks to the extractable velcro separators.
- Reflective strips night and day. - 2 medium-sized external pockets (one of them with an oxygen tube extraction ring)
- 2 large capacity external pockets, with elastics, flat and bulky pockets.
- 5 extractable compartments with the transparent frontal part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: surgical and miscellaneous equipment, yellow: pediatric equipment).
- Auxiliary compartment of gray color.
- Isothermal amp maker with 56 elastics. - Space for tips and masks
- Space for oxygen cylinder max. 5L, with fixing tapes.
- Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana.
- Material: 1000D polyamide
- Color: Red.
- Available in waterproof version
- Weight: 4.12 kg.
- Outdoor measures: 61 x 29.5 x 35 cm / 47.04 L
- Maximum recommended capacity: 15 kg.
- Sold empty

## CRITICAL URGENCY ELITE BAGS, RED

Ref. MECRIT01\_R



### Features

CE Marking	NONE	Medical device	No
Bar code	8436564540171	HS Code	4202921900
Warranty (months)	24	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	31	Bare Unit / Width (cm)	29.5
Packed unit / Height (cm)	24	Bare Unit / Height (cm)	35
Packed unit / Length (cm)	58	Bare Unit / Length (cm)	61
Packed unit / Weight	4	Bare Unit / Weight (kg)	4.12

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---