



HEAD IMMOBILIZER, ORANGE

Ref. SBRAN10_O



11/12/2025 21:40:26 – For information purposes only, can be subject to change

HEAD IMMOBILIZER, ORANGE

Ref. SBRAN10_O



=====

Description

- Manufactured with high density plastic materials, wide holes at the ears to monitor the auditory canal of the patient
- It can be used with a stretcher and a hard board
- Minimum Interference with radiography systems or MRI.
- Plastic Coated, water-resistant and easy to clean
- Dimensions (L x W x H): 40 x 25 x 16 cm
- Weight: 1 kg
- Single colour: orange

HEAD IMMOBILIZER, ORANGE

Ref. SBRAN10_O



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040094656	HS Code	94029000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	26	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	14	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	40
Packed unit / Weight	0.95	Bare Unit / Weight (kg)	1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---