


AC/DC ADAPTOR FOR AIVIA

Ref. SSACCE013



Description

Transformer 24 V DC 75W for cabinet AIVIA

11/12/2025 21:33:09 – For information purposes only, can be subject to change

AC/DC ADAPTOR FOR AVIA

Ref. SSACCE013



Features

CE Marking	NONE	Medical device	No
Bar code	3661041009253	HS Code	901889010
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	6
Packed unit / Length (cm)	10	Packed unit / Weight	0.265
Bare Unit / Weight (kg)	0.245		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--