

Ref. C1.416.114

LUXAMED®



CE



13/12/2025 09:33:42 – For information purposes only, can be subject to change

=====

Ref. C1.416.114

LUXAMED®

=====

Description

Available in several colours:

White (C1.416.914)

Light blue (C1.416.934)

Fuchsia (C1.416.974)

Yellow (C1.416.944)

Kids (C1.816.914)

Black (C1.416.114)

Green (C1.416.954)

Purple (C1.416.964)

Ref. C1.416.114

LUXAMED®

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963417	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.17		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Black
--------	-------