

## ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L

Ref. BCFLOV10

CA-MI



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**IIa**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

## ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L

**CA-MI**

Ref. BCFLOV10

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

### Description

Ca-mi New Askir 20 Phlegm Vacuum - FLOVAC® Single Use Bottle/Bag, 1 L

Fixed electric surgical Suction Unit, designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children.

Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: ENT, Home care, Room, Tracheostomy.

-  
Composed of: New Askir 20 Suction Unit + Flovac 1L Kit (1x Reusable Flovac 1L Jar / 1x Sterile UU pocket 1L Flovac / 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (1 pcs) / 1x Configuration Manual Suction Units

-  
The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Optional roller foot

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

-  
Engine: Oil-free piston pump requiring no maintenance

-  
Power supply: 230V-50Hz - other configurations available on request

- Maximum vacuum (adjustable): 0.75 Bar -75 kPa -563 mmHg

-  
Max free air flow: 16 l/min

-  
Noise level: 60.5 dBA

- Operating

cycle: Uninterrupted

- Weight: 2.5 kg / Dimensions: 35 x 21 x 18 cm

## ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L

Ref. BCFLOV10

**CA-MI**

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040200019	HS Code	90189084
Warranty (months)	24	Country	Italy

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.5
Bare Unit / Weight (kg)	2.5		

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	1	Maximum humidity	85

#### OTHER FEATURES

## ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L

**CA-MI**

Ref. BCFLOV10

### Content

Ref	Designation	Quantity
BCASPI99	DISPOSABLE FLOVAC® LINER, 1000 ML for Aspiret, Askir 20, 30	1
BCASPI164	Male Connector (for Vacuum Port) for Askir 20 pack of 5's	1
BCASPI102	Container Flovac 1L (reusable) for Aspiret, Askir 20 and 30	1
BCASPI504	IFU for CA-MI suction units	1
BCASPI06	New Askir 20 electrical suction unit, 1L autoclavable polycarbonate j	1
BCASPI111	Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30	1