

## ELECTRONIC WRIST BLOOD PRESSURE MONITOR NEOTENS® NOMAD SPENGLER

Ref. NEO100



CE  
0123



15/12/2025 20:05:44 – For information purposes only, can be subject to change

## ELECTRONIC WRIST BLOOD PRESSURE MONITOR NEOTENS® NOMAD SPENGLER



Ref. NEO100

**NEOTENS® NOMAD : compact wrist blood pressure monitor, designed for precision and mobility**

### Description

The **NEOTENS® NOMAD** by Spengler is an electronic wrist blood pressure monitor that combines reliability, comfort, and portability. Its ultra-compact design (70.9 x 70.9 x 26.5 mm) and lightweight build (96g) make it ideal for use in medical offices, home visits, or self-monitoring at home.

Equipped with Spengler's upward measurement technology, it ensures a fast and comfortable reading, while its LCD display guarantees clear visibility. Certified as a Class IIa medical device (CE marking), it meets the highest European safety and performance standards.

- **Display** : clear LCD screen (32.5 x 40.5 mm)
- **Measurement method** : Oscillometric, upward measurement for optimal comfort
- **Measurement site** : Wrist
- **Range** :  
Systolic pressure : 60 - 260 mm Hg  
Diastolic pressure : 40 - 200 mm Hg  
Pulse : 30 - 199 beats/min
- **Accuracy** :  $\pm 3$  mmHg (pressure) ;  $\pm 5$  % (pulse)
- **Memory** : 50 readings per user x 2 users, with a date and time stamp
- **Advanced features** :  
Irregular heartbeat detection  
Low battery indicator  
Automatic shut-off after 1 minute of inactivity
- **Power supply** : 2 x AAA batteries ( $\approx 300$  measurement cycles)
- **Cuff** : Wrist circumference 13.5 - 21.5 cm

# ELECTRONIC WRIST BLOOD PRESSURE MONITOR NEOTENS® NOMAD SPENGLER



## Ref. NEO100

- **Robustness** : IP22 protection rating
- **Weight & Dimensions** : 96 g (without batteries), 70.9 x 70.9 x 26.5 mm
- **Medical classification** : Class IIa device (CE marking)
- **Warranty** : 3 years (excluding consumables)
- **Compatibility** : **SPENGLER HTA mobile app** (manuel entry or photo) for simplified monitoring

This monitor combines technology and ease of use, making it an ideal companion for both healthcare professionals and patients who require accurate and regular monitoring.

### Box contents :

- 1 NEOTENS® NOMAD wrist blood pressure monitor
- 2 x AAA alkaline batteries
- 1 user manual

## Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6942284608388	HS Code	90189010
Warranty (months)	36	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	8	Packed unit / Height (cm)	10
Packed unit / Length (cm)	9	Packed unit / Weight	0.31

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	60
Minimum humidity	NA	Maximum humidity	93

### OTHER FEATURES