


Ref. 337843

 Spengler



CE



 16/12/2025 21:15:48 – For information purposes only, can be subject to change

Ref. 337843



Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Spengler hand-held BPMs Protocol patient monitors

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026943	HS Code	9018901000
Warranty (months)	6	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.125

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---