

## PULSE STETHOSCOPE / RED

Ref. 507522



CE



## PULSE STETHOSCOPE / RED

Ref. 507522



### Description

For daily blood pressure monitoring. Powder-coated or chromium-plated (black model) zinc alloy single head chestpiece. Epoxy resin, high resolution acoustic membrane. Screw-on flexible tips to match the PVC tubing. Angled binaura with built-in spring.

### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 3700446009281 | HS Code                   | 90181990 |
| Warranty (months)           | 36            |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.22     |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Colors                      | Navy blue     | Packed unit / Quantity    | 1        |