


AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS102

DDM




CE



17/12/2025 20:02:12 – For information purposes only, can be subject to change

AUTOCLAVABLE POLYESTER LOWER LIMB CUFF STRAIGHT 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GCS102



Those lower limb cuffs are autoclavable!

Description

Lower limb cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Lower limb cuff infant	35.5 x 6 cm
Lower limb cuff child	53 x 6 cm
Lower limb cuff adult	69 x 9 cm
Lower limb cuff adult L	81 x 9 cm
Lower limb cuff adult XL	85.5 x 12 cm

**AUTOCLAVABLE POLYESTER LOWER LIMB
CUFF STRAIGHT 1 BLADDER FOR DDM
PNEUMATIC TOURNIQUET****Ref. GCS102****Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468203643	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES