

CALL'S ELITE BAGS, BLUE

Ref. MECALL01_B

ELITE BAGS



22/12/2025 20:35:42 – For information purposes only, can be subject to change

CALL'S ELITE BAGS, BLUE

Ref. MECALL01_B



Designed to transport the sanitary equipment of professionals (nurse, chiropodists, veterinarians, doctor ...) during scheduled visits

Description

- Padlocked zipper- 5 extractable compartments with a transparent front side
- Opaque extractable compartment
- Extractable compartment with velcro separators
- Isothermal compartment with 43 elastic bands for medicines- Net pocket
- 4 elastics for fixing the material
- Mini-fastener
- Auxiliary bag for large volumes
- Space for personal items
- Shoulder strap, handle, concealed backpack straps
- Material: nylon 230D
- Color: Blue
- Weight: 2.35 kg
- Outdoor measures: 40 x 30 x 12 cm / 14.4 L
- Maximum recommended capacity: 5 kg
-
- Optional trolley
- Sold empty

CALL'S ELITE BAGS, BLUE

Ref. MECALL01_B



Features

CE Marking	NONE	Medical device	No
Bar code	8436564540058	HS Code	4202121900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	12
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	40
Packed unit / Weight	2.96	Bare Unit / Weight (kg)	2.35

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---