


CRITICAL URGENCY ELITE BAGS, RED WATERPROOF

Ref. MECRIT02_R

ELITE BAGS



22/12/2025 20:35:09 – For information purposes only, can be subject to change

CRITICAL URGENCY ELITE BAGS, RED WATERPROOF



Ref. MECRIT02_R

Description

Allows all equipment intended for advanced first aid of emergency services to be transported, facilitating access to equipment quickly and efficiently.

Interior totally adaptable to the needs of the user thanks to the extractable velcro dividers.

Reflective strips night and day.

2 medium-sized external pockets (one of them with an oxygen tube extraction ring)

2 large capacity external pockets, with elastics, flat and bulky pockets.

5 extractable compartments with transparent frontal part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: surgical and miscellaneous equipment, yellow: pediatric equipment).

Auxiliary compartment of gray color.

Isothermal amp maker with 56 elastics. Space for tips and masks

Space for oxygen cylinder with fixing tapes.

Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana.

Colour: Red.

Weight: 4.12 kg.

Outdoor measures: 61 x 29.5 x 35 cm / 47.04

I. Maximum recommended capacity: 15 kg.

CRITICAL URGENCY ELITE BAGS, RED WATERPROOF



Ref. [MECRIT02_R](#)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040167473	HS Code	4202921900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	29.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	35
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	61
Packed unit / Weight	nc	Bare Unit / Weight (kg)	4.12

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---