

BAG FOR OXYGEN CYLINDER FROM 1 TO 3L, RED EB02.037

Ref. METUBE3L_R

ELITE BAGS



22/12/2025 20:31:44 – For information purposes only, can be subject to change

BAG FOR OXYGEN CYLINDER FROM 1 TO 3L, RED EB02.037



Ref. METUBE3L_R

Description

- Cylindrical bag for oxygen therapy
- Sold empty
- For oxygen bottle of 1 to 3 liters
- Measures: 46 x 20 x 15 Ø cm; weight 0.72 Kg
- Capacity: 10 L
- Material: 1000D Polyester. Washable.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040202921	HS Code	4202921900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Bare Unit / Width (cm)	15
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	46	Bare Unit / Length (cm)	46
Packed unit / Weight	0.72	Bare Unit / Weight (kg)	0.72

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---