

=====

Ref. E1.421.012

**LUXAMED®**



CE



=====

## Description

=====

 25/12/2025 19:49:58 – For information purposes only, can be subject to change

Ref. E1.421.012

LUXAMED®

## Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 4260381960584 | HS Code                   | 90189084 |
| Country                     | Germany       |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | nc       |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | 5             | Max temperature           | 40       |
| Minimum humidity            | 10            | Maximum humidity          | 90       |
| OTHER FEATURES              |               |                           |          |
| Size                        | Size 1        |                           |          |