

Ref. 522713

 Spengler



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
5 ANS
5
5 YEARS
GUARANTEE



 29/12/2025 22:04:09 – For information purposes only, can be subject to change

Ref. 522713



Description

Cuff available in :

Adult (M), nylon, velcro, black (522711)

Adult (L), nylon, black (522712)

Child/small adult (S), nylon, black (522713)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028190	HS Code	9018199560
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.071
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		