


SMARTLIGHT OTOSCOPE BLACK

Ref. 570315



05/01/2026 22:20:29 – For information purposes only, can be subject to change

SMARTLIGHT OTOSCOPE BLACK

Ref. 570315



The otoscope with conventional lighting

Description

Conventional lighting. Lightweight and robust, this otoscope is perfect for ENT auscultations especially on the move. The XHL lamp diffuses a concentrated and homogenous light beam (optional LED lamp available). Fit with an ABS and Brass head, SMARTLIGHT also comes with end clamp and built-in switch, LR6 AA batteries included.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041557	HS Code	9018908400
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Packed unit / Height (cm)	3
Packed unit / Length (cm)	10	Packed unit / Weight	0.23
Bare Unit / Weight (kg)	0.1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--