

Ref. C1.416.914

LUXAMED®



CE



05/01/2026 22:47:47 – For information purposes only, can be subject to change

Ref. **C1.416.914**

**LUXAMED®**

## Description

Available in several colours:

White (C1.416.914)

Light blue (C1.416.934)

Fuchsia (C1.416.974)

Yellow (C1.416.944)

Kids (C1.816.914)

Black (C1.416.114)

Green (C1.416.954)

Purple (C1.416.964)

## Features

| CE Marking        | CE             | Medical device | Class I  |
|-------------------|----------------|----------------|----------|
| Bar code          | 04260381963868 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Germany  |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | 5  | Max temperature  | 40 |
| Minimum humidity    | 10 | Maximum humidity | 90 |

### OTHER FEATURES

|        |       |
|--------|-------|
| Colors | White |
|--------|-------|

05/01/2026 22:47:47 – For information purposes only, can be subject to change