

RESCUE BACKPACK EMERGENCY, BLUE EM07.004

Ref. MEEMER05_R

EMERGENCY'S™



05/01/2026 23:27:47 – For information purposes only, can be subject to change

RESCUE BACKPACK EMERGENCY, BLUE EM07.004

Ref. MEEMER05_R

EMERGENCY'S™

Its foldable flap opening provides access to the tools necessary for emergency actions

Description

- Main ergonomic handle with Velcro- Large pockets on both sides
- Padded shoulder straps and backs, adjustable backpack
- Large front pocket, ready to accommodate a defibrillator
- Spacious interior for convenient and easy access to the material
- Internal compartments with mesh and transparent windows for the organization of the material
- External dimensions: 37 x 45 x 21 cm / 30L
- Material: 600D Polyester and PVC
- Weight: 1.22 kg
- Available in waterproof version
- Sold empty

Available in :

Blue : MEEMER05_B

Red : MEEMER05_R

RESCUE BACKPACK EMERGENCY, BLUE EM07.004

Ref. [MEEMER05_R](#)

EMERGENCY'S™

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540829	HS Code	420299000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	31	Bare Unit / Width (cm)	37
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	21
Packed unit / Length (cm)	48	Bare Unit / Length (cm)	45
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.22

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Red	Packed unit / Quantity	1
--------	-----	------------------------	---