

ELITE BAGS - PRACTI - BLUE

Ref. MEPRACT02_N

ELITE BAGS



05/01/2026 23:33:25 – For information purposes only, can be subject to change

ELITE BAGS - PRACTI - BLUE

Ref. MEPRACT02_N



Health assistance bag for home consultations or emergency.

Description

Its excellent internal distribution makes it possible to immediately identify materials.

-

2 extractable compartments with transparent front part

-

Isothermal amp maker

-

External dividers

-

Inner and outer pockets- Measures: 40 x

30 x 16 cm

-

Capacity: 19.20 L

-

Weight: 1.82 kg

-

Material: 420D polyester

-

Washable

-

Maximum recommended capacity: 4 kg

- Sold empty

ELITE BAGS - PRACTI - BLUE

Ref. MEPRACT02_N



Features

CE Marking	NONE	Medical device	No
Bar code	8436564541819	HS Code	42029900
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	31	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	43	Bare Unit / Length (cm)	40
Packed unit / Weight	2	Bare Unit / Weight (kg)	1.82
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black		