

## STETHOSCOPE PULSE II SINGLE

Ref. 507701



CE



## STETHOSCOPE PULSE II SINGLE

Ref. 507701



The essential stethoscope. Minimum weight, maximum performance.

### Description

#### HIGH ACOUSTIC QUALITY

- Single tube with high-performance conduit
- High resolution binaural
- High resolution diaphragm

#### COMFORT

- Anti-cold ring
- Silicone ear tips
- New headphone design more adapted to the anatomy
- Aluminium earhook to reduce weight

Available in Carbon (507702), Blueberry (507700), Coral (507701) or Powder Pink (507703).

Available with double bell (adult bell and bell).

## STETHOSCOPE PULSE II SINGLE

Ref. 507701



### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | Class I  |
| Bar code                    | 3700446029876 | HS Code                   | 90181990 |
| Warranty (months)           | 60            |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.226    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | NA            | Max temperature           | NA       |
| Minimum humidity            | NA            | Maximum humidity          | NA       |
| OTHER FEATURES              |               |                           |          |
| Colors                      | Coral         | Packed unit / Quantity    | 1        |