

TROLLEY FOR 5L GAS BOTTLES

Ref. 880030



TROLLEY FOR 5L GAS BOTTLES

Ref. 880030



Description

Trolley for medical gas cylinders.
This trolley makes it possible to transport RDI bottles with confidence.
Compatible 5L MEOPA and O2 RDI bottles

- Dimensions: L300mm*W270mm*H(890-1030)mm
- Weight: 1.764KG/per set

Packaging Details:

Features

CE Marking	NONE	Medical device	No
Bar code	3661040215914	HS Code	90189080
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	60	Packed unit / Height (cm)	31
Packed unit / Length (cm)	31	Packed unit / Weight	2.15
Bare Unit / Weight (kg)	1.764		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		