

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)

Ref. 519050



CE
0459



12/01/2026 21:35:52 – For information purposes only, can be subject to change

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519050

=====

Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

CUFF

- Cotton CUFF, velcro fastening
- Size XXS (10-15 cm), XS (15-22 cm), S (20-28 cm, ref : 519029), M (26-34 cm, ref 519012), or L (32-42 cm), S-M-L (ref :519050)
- Available in grey, navy, red or patterns
- Supplied with carrying case

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)

Ref. 519050



Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024789	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	8
Packed unit / Length (cm)	26	Packed unit / Weight	1.15
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		