

## LITESTICK PENLIGHT SPENGLER

Ref. 570620



## LITESTICK PENLIGHT SPENGLER

Ref. 570620



Pen-shaped medical light, designed for close-up clinical examinations.

### Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009830	HS Code	9018908400
Warranty (months)	24		

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.015		

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	40
---------------------	----	-----------------	----

#### OTHER FEATURES

Colors	Black	Packed unit / Quantity	1
--------	-------	------------------------	---