


## LITESTICK PENLIGHT SPENGLER

Ref. 570623



## LITESTICK PENLIGHT SPENGLER

Ref. 570623



Pen-shaped medical light, designed for close-up clinical examinations.

### Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

### Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | CE            | Medical device            | Class I    |
| Bar code                    | 3700446009847 | HS Code                   | 9018908400 |
| Warranty (months)           | 24            |                           |            |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | 5             | Packed unit / Height (cm) | 5          |
| Packed unit / Length (cm)   | 20            | Packed unit / Weight      | 0.046      |
| Bare Unit / Weight (kg)     | 0.015         |                           |            |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | 10            | Max temperature           | 40         |
| OTHER FEATURES              |               |                           |            |
| Colors                      | Red           | Packed unit / Quantity    | 1          |